

How to Use the Data Exchange

Only for **pathway owners** using a **pro one™** or **Enterprise** system.

You will need to be set up as a user on both the source and destination sites with identical custom lists. These steps will need to be done only by consultants, delegates will not be able to do this.

Click the Settings cog on your pro one or Enterprise system.

Patient	Pathway	C	P
DEMO1, Demo	Knee Pathway (Right)	3	0
DEMO2, Demo	Shoulder Pathway (Left)	1	5
DEMO3, Demo	T + O Baseline Generic Pathway (Unknown)	1	0

Click Sync.

Pathway Type	Status	Action
Generic Pathway	ON	SET DEFAULT
Patient Experience Pathway	OFF	
Lower Limb (Distal Tibia/Foot & Ankle) Tumour/Salvage/Reconstruction Pathway	OFF	
Lower Limb (Femur/Tibia) Tumour/Salvage/Reconstruction Pathway	OFF	
Shoulder And Upper Arm Tumour/Salvage/Reconstruction Pathway	OFF	
Ankle Fracture Pathway (Trauma)	OFF	
Foot and Ankle (Achilles Rupture) Pathway (Trauma)	OFF	

Click Authenticate for the site you would like to sync to. If you would like to sync to more than one site, please repeat this process for each site.

Gold CONSULTANT, DEMO ▾

User Settings 🏠 📄 + 🔍 📊 📄 7 ⚙️ ?

PATHWAY TYPES
9 ITEMS SELECTED

CUSTOM LISTS DELEGATES PROFILE SECURITY SYSTEM SETTINGS NATIONAL IDS SYNC

SYNC JOBS

SYNCABLE PATHWAY LIST SYNC HISTORY

Foot and Ankle Pathway

Destination	Last Run Date
Authenticate Amplitude Demo SA	N/A
Authenticate ClinicalRegistryTest	N/A

Generic Pathway

Destination	Last Run Date
Authenticate Amplitude Demo SA	N/A
Authenticate ClinicalRegistryTest	N/A

T + O Baseline Generic Pathway

Destination	Last Run Date
Authenticate Amplitude Demo SA	N/A
Authenticate MultiTenantTest	N/A

You will then be asked to agree to the terms of the data exchange.

User Settings 🏠 📄 + 🔍 📊 📄 7 ⚙️ 🔧 ?

PATHWAY TYPES
6 ITEMS SELECTED

CUSTOM LISTS DELEGATES PROFILE SECURITY

BOFAS

By setting up this data transfer you are agreeing to the patient data relating to your patients, on designated pathways, being transferred to the nominated clinical registry/ies and/or your individual pro one™ account. As the Data Controller, you confirm you have obtained all the necessary patient consents and abide by the Amplitude Data Processing Agreement. This data transfer will remain active until it is cancelled by you.

Cancel **Agree & Continue**

Enter your login credentials for the destination site.

Click Log in.

User Settings 🏠 📄 + 🔍 📊 📄 7 ⚙️ ?

PATHWAY TYPES
10 ITEMS SELECTED

CUSTOM LISTS DELEGATES PROFILE SECURITY SYSTEM SETTINGS NATIONAL IDS SYNC

CLINICALREGISTRYTEST

Please login to the destination system:

Amplitude Development Environment

GOLD (TEST)
For Internal Use Only

Log in

Username
DemoConsultant

Password
.....

Best friend at school
|

Cancel **Log in**

Switch on the pathways you would like to sync to that site.

The scheduler will now sync these pathways at the allotted time.

User Settings

SYNC JOBS

SYNCABLE PATHWAY LIST | SYNC HISTORY

Knee (Primary ACL Reconstruction) Pathway (NLR)

Destination	Last Run Date
ON NLR	11 Sep 2019 06:01

Knee (Primary ACL Repair) Pathway (NLR)

Destination	Last Run Date
ON NLR	11 Sep 2019 06:01

Knee (UKKOR - Osteotomy) Pathway

Destination	Last Run Date
ON UKKOR	11 Sep 2019 06:01

To check your data transfers click Sync History. Amplitude recommends checking the list regularly to make sure all data sync's and patients are sending successfully.

User Settings

SYNC JOBS

SYNCABLE PATHWAY LIST | SYNC HISTORY

Show All Destinations | Show All Pathway Types

Destination	Pathway Type	Start Date	Status	✓	!	Show Records
NLR	Knee (Primary ACL Reconstruction) Pathway (NLR)	11 Sep 2019 06:01	🔄	37	3	Show Records
NLR	Knee (Primary ACL Repair) Pathway (NLR)	11 Sep 2019 06:01	🔄	4	1	Show Records
UKKOR	Knee (UKKOR - Osteotomy) Pathway	11 Sep 2019 06:00	🔄	59	5	Show Records
NLR	Knee (Primary ACL Reconstruction) Pathway (NLR)	04 Sep 2019 06:01	🔄	36	3	Show Records
NLR	Knee (Primary ACL Repair) Pathway (NLR)	04 Sep 2019 06:01	🔄	4	1	Show Records
UKKOR	Knee (UKKOR - Osteotomy) Pathway	04 Sep 2019 06:00	🔄	57	5	Show Records
NLR	Knee (Primary ACL Reconstruction) Pathway (NLR)	28 Aug 2019 06:01	🔄	36	3	Show Records
NLR	Knee (Primary ACL Repair) Pathway (NLR)	28 Aug 2019 06:01	🔄	4	1	Show Records
UKKOR	Knee (UKKOR - Osteotomy) Pathway	28 Aug 2019 06:00	🔄	57	5	Show Records
NLR	Knee (Primary ACL Reconstruction) Pathway (NLR)	21 Aug 2019 06:00	🔄	33	3	Show Records
NLR	Knee (Primary ACL Repair) Pathway (NLR)	21 Aug 2019 06:00	🔄	4	1	Show Records

You can choose the site and pathway to see the history of the sync job. Patients that fail to send appear in the column. Click 'Show Records' to see who and why they failed to send. If you see an error with the data exchange please contact Amplitude Customer Support.